

Fecal Occult Blood Collection Instructions

1650 Cowles (Fairbanks Imaging Center) Fairbanks, AK 99701 907-458-5640 Or

1101 and 1001 Noble (in the TVC building) Fairbanks, AK 99701 907-459-3500, ext. 2229

Drug Guidelines:

- For seven days before and during the stool collection period avoid non-steroidal antiinflammatory drugs such as ibuprofen, naproxen or aspirin (more than one adult aspirin a day).
- Acetaminophen can be taken as needed.
- For **three** days before and during stool collection period avoid vitamin C in excess of 250 mg a day from supplements and citrus fruits and juices. 100% of RDA of vitamin C for an adult is 60 mg a day. Some iron supplements contain vitamin C in excess of 250 mg.

Note: Please talk to your physician or pharmacist if you have any questions about mediations you take regularly.

Diet Guidelines

- For three days before and during stool collection period avoid red meats (beef, lamb and liver).
- Eat a well balanced diet including fiber such as bran cereals, fruit and vegetables.

General Guidelines:

- Do not collect sample if blood is visible in your stool. **Contact your doctor**.
- For the most accurate test results collect each stool sample before contact with water. You may use any clean, dry container for collection.

Collection:

- 1. Remove slide. Using ball point pen, write your name (first and last) and date of birth on the front of the slide.
- 2. Fill in sample collection date on section before a bowel movement.
- 3. You may use any clean, dry container to collect your sample. Let stool fall into container.
- 4. Use one stick to collect a small sample. Apply a thin smear covering box A.
- 5. Collect a second sample from a different part of the stool sample with the same stick. Apply a thin smear covering box B.
- 6. Discard stick in a waste container.
- 7. Close and secure front flap section by inserting it under tab. Allow sample to dry.

Transportation:

Return completed slide no later than 14 days after collection.

If you have any questions please call Fairbanks Memorial Hospital Laboratory at 907-458-5640 or your practitioner.

Fecal Occult Blood Collection Instructions